

Zavation Interbody Fusion Devices

Cervical




Heights: 5mm-12mm
 Widths: 14mm
 Lengths: 12mm
 Lordotic: 6°, 10°

Lateral Lumbar




8mm-16mm
 19mm, 23mm
 45mm-65mm
 Parallel, 7°

Anterior Lumbar



9mm-21mm
 30mm 36mm
 24mm 28mm
 8°, 15°

Posterior Lumbar



TLIF

Heights: 7mm-16mm
 Widths: 10mm
 Lengths: 25 and 30mm




Parallel
 T-PLIF
 PLIF(23mm)

7mm-16mm
 9mm, 11mm
 23mm, 27mm, 32mm



Lordotic
 T-PLIF-L
 PLIF-L(23mm)

8mm-16mm
 9mm, 11mm
 23mm, 27mm, 32mm



Convex
 T- PLIF-C
 PLIF-C(23mm)

8mm-16mm
 9mm, 11mm
 23mm, 27mm, 32mm



Rotatable
 T-PLIF-R
 PLIF-R(23mm)

10mm-16mm
 9mm, 11mm
 23,mm,27mm, 32mm

*23mm devices are to be used bi-laterally

Benefits and Features

- PEEK cages with large graft hole and robust insertion tool
- Threaded central hole with anti-rotation holes for inserter
- Bulleted distal ends for easier insertion
- Tantalum marker



LNE G
 MED
 ISO 13485:2003
 Certificate No: 22592rev 0

Brandon, MS 39047
 501 Avalon Way
 Ph: 601-919-1119
www.zavation.com